


Flasks holders and supports/volumetric

For product assistance please telephone: 01509 231166



Flask supports, closed lead rings

02



Closed lead rings coated with chemical resistant Vikem®. Ideal for stabilising Erlenmeyer, volumetric and straight necked flasks.

Catalogue No	Alt. No	Mass, g	Diameter (internal), mm	Diameter (external), mm	For flask capacity, mL
FHM-415-010F	F18307-0003	165	19	51	5 to 10 cylinders
FHM-415-030W	F18307-0005	237	43	76	125 to 500 flasks
FHM-415-050Q	F18307-0010	454	51	89	250 to 1,000 flasks
FHM-415-070K	F18307-0015	681	56	102	500 to 2,000 flasks
FHM-415-090E	F18307-0020	908	71	114	1,000 to 4,000 flasks



Flask supports, open lead rings

02



Open lead rings coated with chemically resistant Vikem®. Ideal for stabilising Erlenmeyer, volumetric and straight necked flasks.

Catalogue No	Alt. No	Mass, g	Diameter (internal), mm	Flask size, mL
FHM-410-010M	F18308-0500	230	43	125 to 500
FHM-410-030G	F18308-1000	450	51	250 to 1,000
FHM-410-050A	F18308-2000	680	56	500 to 2,000
FHM-410-070R	F18308-4000	900	71	1,000 to 4,000



Flasks, borosilicate, volumetric, class A, small capacity

01



NEW

Mini volumetric flasks to facilitate handling of small samples. Manufactured in chemically resistant Pyrex® glass, with white enamel graduations and polyethylene stoppers.

Catalogue No	Alt. No	Capacity, mL	Tolerance, mL	Stopper size	Pack qty
FHM-395-010D	3151/01	1	0.02	7/16	5
FHM-395-020A	3151/02	2	0.02	7/16	5
FHM-395-030U	3151/04	5	0.02	10/19	5



Flasks, borosilicate glass, volumetric, trapezoidal, Class A

01



Mini volumetric flasks with tight accuracy tolerances, and shaped to maximise stability while allowing access to contents. Individually calibrated and tested.

Catalogue No	Alt. No	Capacity, mL	Tolerance, mL	Pack qty
Clear				
FHM-440-010X	RB257/A/5	1	0.010	5
FHM-440-030Y	RB257/B/5	2	0.015	5
FHM-440-050C	RB257/C/5	3	0.025	5
FHM-440-070T	RB257/D/5	4	0.025	5
FHM-440-090N	RB257/E/5	5	0.025	5
Amber				
FHM-450-010H	RB257/AM/A/5	1	0.010	5
FHM-450-030B	RB257/AM/B/5	2	0.015	5
FHM-450-050S	RB257/AM/C/5	3	0.025	5
FHM-450-070M	RB257/AM/D/5	4	0.025	5
FHM-450-090G	RB257/AM/E/5	5	0.025	5